# Beißring SU

DE

EN

FR

IT

ES

NL

CZ

PT

SE

NO

DK

RU

CN

KR

TU

PL

SA

EL

HR

RO

<b>DE</b>	Beißring	Gebrauchsanweisung	<b>3</b>
<b>EN</b>	Bite Block	Instruction Manual	<b>8</b>
<b>FR</b>	Ouvre-bouche	Manuel d'utilisation	<b>13</b>
<b>IT</b>	Boccaglio	Manuale di istruzioni	<b>18</b>
<b>ES</b>	Anillo de dentición	Manual de instrucciones	<b>23</b>
<b>NL</b>	Bijtvergrendelin	Handleiding	<b>28</b>
<b>CZ</b>	Kousátka	Návod k použití	<b>33</b>
<b>PT</b>	Bloqueio de mordidas	Manual de Instruções	<b>38</b>
<b>SE</b>	Bittlås	Instruktionsbok	<b>43</b>
<b>NO</b>	Bittlås	Instruksjonshåndbok	<b>48</b>
<b>DK</b>	Bitlås	Betjeningsvejledning	<b>53</b>
<b>RU</b>	Замок для укусов	Руководство по эксплуатации	<b>58</b>
<b>CN</b>	出牙环	说明书	<b>63</b>
<b>KR</b>	물린 잠금	사용 설명서	<b>67</b>
<b>TU</b>	Diş kaşığıcı	Kullanım talimatları	<b>72</b>
<b>PL</b>	Pierścienie uszczelniające	Instrukcje użytkowania	<b>77</b>
<b>SA</b>	قضاضع	إدخيتس الل تامي لع	<b>82</b>
<b>EL</b>	Μανταλάκι δαγκώματος	Οδηγίες χρήσης	<b>87</b>
<b>HR</b>	Zaključavanje ugriza	Uputa za uporabu	<b>92</b>
<b>RO</b>	Inele de dentiție	Instrucțiuni de utilizare	<b>97</b>

## Symbole

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## Symbols

EN

	Medical Device		Latex-Free
	Catalog Number		non-sterile
	Batch		
	Warning! Read the Instructions for use before use		
	Usable until		
	For single use only. Do not reuse		
	Do not re-sterilize!		
	Keep dry		
	Protect from sunlight		
	Manufacturer		
	Date Of Manufacture		
	Do not use if packaging is damaged		

---

## Copyright

The copyright to these instructions for use remains with the company ENDO-FLEX GmbH. All rights, including those of translation, are reserved. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), in whole or in part, without the written permission of the company ENDO-FLEX GmbH or processed, reproduced, translated or distributed using electronic systems. Violations can result in criminal penalties.

## Limitation Of Liability

The does ENDO-FLEX GmbH not assume any liability for personal injury, damage to property, damage to the product or consequential damage caused by non-compliance with these instructions for use, improper use of the product or any other action by unqualified and not ENDO-FLEX GmbH trained personnel. It is also forbidden to carry out modifications or technical modifications to the product by yourself.

## Purpose of the Instructions for use

These Instructions for use are intended for use by the medical personnel working with the product. It describes the structure, function and operation of the product. Read these Instructions for use carefully and make sure that you understand them and that they are accessible to healthcare professionals do it!

Follow the instructions to ensure safe working. Also, follow the instructions in other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The use of the product may only be carried out by qualified and, if necessary, on request ENDO-FLEX GmbH by trained medical personnel. The use of this product by healthcare professionals is always carried out taking into account the individual benefits of the patient at their risk.

Always have a copy of these Instructions for use at hand when using.

## Classification

The Beißring is designated as a class I product according to Medical Device Regulation (MDR) – (2017/745) via medical devices.

## Reporting Obligation


The indicates ENDO-FLEX GmbH that the user is obliged to report to the competent authority of his country and at the same time all adverse events occurring during the use of this product and concerning this ENDO-FLEX GmbH product.

## Language

These and all other instructions for use are available as a PDF download at [www.endo-flex.de](http://www.endo-flex.de) .

## 1. Product Description

The ENDO-Flex GmbH offers teethers with fixation tape in different versions for adults as well as for children. These include teethers in standard and special versions with tongue plates.

This Instructions for use is valid for products listed below:

- 600005
- 600006
- 600010
- 600020

### 1.1. Contents of the packaging

600005 / 600006:

- 1 Carton
- 5 Beißring SU (single use), individually non-sterile packed
- 1 Instructions For Use

600010 / 600020:

- 1 outer carton (dispenser box of 100 pieces)
- 100 Beißring SU (single use), individually non-sterile packed
- 1 Instructions For Use

### 1.2. Intended Use

Beißring Are used to provide a safe passage for endoscopes and/or esophageal dilators, in which they protect them from damage/destruction by unintended upsetting of the patient to be treated. Teethers SU are equipped with a head strap for secure fixation.

#### 1.2.1 Duration Of Use

The application duration is specified for less than 60 minutes.

#### 1.2.2 Lifetime

The life of the product is 3 years.

#### 1.2.3 Indication

- Safe passage of endoscopes/esophageal dilators

#### 1.2.4 Contraindications

- Not known

#### 1.2.5 Potential Complications/Risks

- Not known

#### 1.2.6 Target group of users

Trained specialists and specialists.

### 1.3. Sterility

These disposable medical devices are supplied non-sterile. Treatment and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

## 2. Safety Information

The safety and warning information used in these Instructions for use is based on ANSI Z535 standards.

<b>DANGER</b>	Warns of an imminent danger that leads to death or serious injury.
<b>WARNING</b>	Warns of a potentially dangerous situation that could lead to death or serious injury.
<b>CAUTION</b>	Warns of a potentially dangerous situation that could lead to moderate or minor personal injury.

<b>NOTE</b>	Warns of a potentially dangerous situation that could lead to property damage and environmental damage.
<b>SAFETY INSTRUCTIONS</b>	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

- Undamaged packaging.
- Damage to the product (cracks on the tube, kinks, deformations).
- Check products for correct surfaces, correct assembly and functionality.
- Products that have not passed the functional test must not be used as their sterility and product safety are not guaranteed. Dispose of the products to be complained about accordingly or return them to the manufacturer.

## 2.1. General

Follow these instructions, as well as instructions from the compatible components and hospital-internal regulations for infection prevention, safe use, cleaning and sterilization.

- All components should be carefully checked for compatibility and integrity prior to use. Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.
- Never use the product outside the recommended technical specifications (intended use).
- Never manipulate the structural conditions of the instrument, avoid kinks and other damage, immediately discontinue use in the event of a malfunction!
- It is absolutely necessary to wear protective clothing (gloves, mouthguards, goggles, gown, etc.)!

<b>NOTE</b>
Never use the product outside the recommended technical specifications.

 <b>WARNING</b>
Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.

<b>SAFETY INSTRUCTIONS</b>
Take appropriate measures to protect the patient and yourself from contamination and injury. Follow national and internal hygiene and safety guidelines. If in doubt, contact your local representative.

 <b>WARNING</b>
Never manipulate the structural conditions of the instrument; in the event of a malfunction, discontinue use.

## 3. Application

### 3.1. Preparation

The medical devices must be tested prior to use with regard to the following aspects:

- Expiration Date.

### 3.1.1 Work Preparation & Testing

<b>NOTE</b>
Do not use products that have expired.

<b>NOTE</b>
If the packaging is damaged before use, replace the instrument with a new one.

<b>NOTE</b>
If you notice any irregularities, replace the instrument with a new one.

 <b>WARNING</b>
Using a defective instrument can cause damage to the endoscope and result in injury to patients or users.

### 3.2. Application

1. Unpack the teether.
2. Position the teether between the patient's upper and lower jaw.

<b>NOTE</b>
Take care not to damage the teeth and tongue!

3. Place the fixation strap attached to the teether on one side around the patient's head and attach it to the opposite side of the teether by inserting the side fixation strap bracket through one of the holes in the tape.

<b>NOTE</b>
Make sure that the fuser belt is not too tight!

### 3.3. Preparation

#### WARNING

Reconditioning and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

### 4. Technical Data

The list of all variants as well as the materials and individual parts used can be found in the technical product data sheet and can be obtained on request from ENDO-FLEX GmbH the and/or on the website of the . ENDO-FLEX GmbH

### 5. Storage & Transportation

This product may only be transported and stored in the packaging provided for it. In addition, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store instruments near chemicals, disinfectants, or radioactive radiation.

#### SAFETY INSTRUCTIONS

Do not place objects on the storage packaging and the sterile barrier system!

### 6. Disposal

After use, this product may present a biological hazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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