
Endoskop-Reinigungsbürsten

DE

EN

FR

IT

ES

NL

CZ

PT

SE

NO

DK

RU

CN

KR

TU

PL

SA


EL

HR

RO

DE	Endoskop-Reinigungsbürste	Gebrauchsanweisung	3
EN	Endoscope Cleaning Brush	Instruction Manual	8
FR	Brosse de nettoyage d'endoscope	Manuel d'utilisation	13
IT	Spazzola per la pulizia dell'endoscopio	Manuale di istruzioni	18
ES	Cepillo de limpieza para endoscopios	Manual de instrucciones	23
NL	Endoscoop reinigingsborstel	Handleiding	28
CZ	Čisticí kartáč na endoskop	Návod k použití	33
PT	Escova de limpeza de endoscópio	Manual de Instruções	38
SE	Endoskop rengöringsborste	Instruktionsbok	43
NO	Endoskoprensebørste	Instruksjonshåndbok	48
DK	Endoskoprensebørste	Betjeningsvejledning	53
RU	Щетка для очистки эндоскопа	Руководство по эксплуатации	58
CN	内窥镜清洁刷	说明书	63
KR	내시경 청소용 브러시	사용 설명서	68
TU	Endoskop temizleme fırçası	Kullanım talimatları	73
PL	Szczotka do czyszczenia endoskopu	Instrukcje użytkowania	78
SA	راظنملا فبيظنت ءاشرف	مءدخستالل تاميلاءع	83
EL	Βούρτσα καθαρισμού ενδοσκοπίου	Οδηγίες χρήσης	88
HR	Endoskop četka za čišćenje	Uputa za uporabu	93
RO	Perii de curățare a endoscopului	Instrucțiuni de utilizare	98

Symbole

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Symbols

EN

	Medical Device		Working Length
	Catalog Number		Outer Diameter
	Batch		Minimum Working Channel
	Warning! Read the Instructions for use before use		non-sterile
	Usable until		
	For single use only. Do not reuse		
	Do not re-sterilize!		
	Keep dry		
	Protect from sunlight		
	Manufacturer		
	Date Of Manufacture		
	Do not use if packaging is damaged		

Copyright

The copyright to these instructions for use remains with the company ENDO-FLEX GmbH. All rights, including those of translation, are reserved. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), in whole or in part, without the written permission of the company ENDO-FLEX GmbH or processed, reproduced, translated or distributed using electronic systems. Violations can result in criminal penalties.

Limitation Of Liability

The does ENDO-FLEX GmbH not assume any liability for personal injury, damage to property, damage to the product or consequential damage caused by non-compliance with these instructions for use, improper use of the product or any other action by unqualified and not ENDO-FLEX GmbH trained personnel. It is also forbidden to carry out modifications or technical modifications to the product by yourself.

Purpose of the Instructions for use

These Instructions for use are intended for use by the medical personnel working with the product. It describes the structure, function and operation of the product. Read these Instructions for use carefully and make sure that you understand them and that they are accessible to healthcare professionals do it!

Follow the instructions to ensure safe working. Also, follow the instructions in other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The use of the product may only be carried out by qualified and, if necessary, on request ENDO-FLEX GmbH by trained medical personnel. The use of this product by healthcare professionals is always carried out taking into account the individual benefits of the patient at their risk.

Always have a copy of these Instructions for use at hand when using.

Classification

The Endoscope Cleaning Brush is designated as a class I product according to Medical Device Regulation (MDR) – (2017/745) via medical devices.

Reporting Obligation

The indicates ENDO-FLEX GmbH that the user is obliged to report to the competent authority of his country and at the same time all adverse events occurring during the use of this product and concerning this ENDO-FLEX GmbH product.

Language


These and all other instructions for use are available as a PDF download at www.endo-flex.de .

1. Product Description

The ENDO-Flex GmbH offers endoscope-Reinig in different versions:

Double-sided endoscope cleaning brushes for manual cleaning of the inner sides of the working or suction channel of commercially available endoscopes. Cleaning brushes are available in different lengths and diameters adapted to the endoscopes. The bristles and shaft of the cleaning brushes are made of high-quality plastics.

Double-sided valve cleaning brushes for the manual cleaning of the inside of the proximal biopsy, suction and flushing valve on the endoscope as well as the distal albarrane lever. The bristles and shaft of the cleaning brushes are made of high-quality plastics.

This Instructions for use is valid for the products/product groups listed below:

For Work Channels:

- EBS41181-A.
- EBS41221-C.
- EBS41221-G.
- NEBS41221-C.
- NEBS41221-G.
- NEBS41221-M.

For Valves:

- ZB0050

In the kit (for working channels and valves):

- NEBS41221-C-SET01
- NEBS41221-G-SET01

1.1. Contents of the packaging

For working channels - EBS series:

- 1 Carton
- 10 Cleaning brushes SU (SingleUse) non-sterile packed
- 1 Instructions For Use

For working channels - NEBS series:

- 1 Carton
- 25 Cleaning brushes SU (SingleUse) non-sterile packed
- 1 Instructions For Use

For valves - ZB series:

- 1 Carton
- 25 Cleaning brushes SU (SingleUse) non-sterile packed
- 1 Instructions For Use

In Set:

- 1 Carton
- 25 Cleaning brush sets SU (SingleUse) non-sterile packed
- 1 Instructions For Use

1.2. Intended Use

Endoscope cleaning brushes are used for manual pre-cleaning of work/suction channels and the suction, flushing and biopsy valve immediately after use of the endoscope to prevent tissue residues from drying in. They are part of a hospital/practice-internal treatment procedure.

1.2.1 Duration Of Use

The application duration is specified for less than 60 minutes.

1.2.2 Lifetime

The life of the product is 3 years.

1.2.3 Indication

- Manual cleaning of contaminated endoscope channels
- Manual cleaning of contaminated endoscope valves

1.2.4 Contraindications

- Not known

1.2.5 Potential Complications/Risks

- Not known

1.2.6 Target group of users

Trained specialists and specialists.

1.3. Sterility

These disposable medical devices are supplied non-sterile. Treatment and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

2. Safety Information

The safety and warning information used in these Instructions for use is based on ANSI Z535 standards.

DANGER	Warns of an imminent danger that leads to death or serious injury.
WARNING	Warns of a potentially dangerous situation that could lead to death or serious injury.
CAUTION	Warns of a potentially dangerous situation that could lead to moderate or minor personal injury.
NOTE	Warns of a potentially dangerous situation that could lead to property damage and environmental damage.
SAFETY INSTRUCTIONS	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

2.1. General

Follow these instructions, as well as instructions from the compatible components and hospital-internal regulations for infection prevention, safe use, cleaning and sterilization.

- All components should be carefully checked for compatibility and integrity prior to use. Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.
- Never use the product outside the recommended technical specifications (intended use).
- Never manipulate the structural conditions of the instrument, avoid kinks and other damage, immediately discontinue use in the event of a malfunction!
- It is absolutely necessary to wear protective clothing (gloves, mouthguards, goggles, gown, etc.)!

NOTE

Never use the product outside the recommended technical specifications.

WARNING

Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.

SAFETY INSTRUCTIONS

Take appropriate measures to protect the patient and yourself from contamination and injury. Follow national and internal hygiene and safety guidelines. If in doubt, contact your local representative.

WARNING

Never manipulate the structural conditions of the instrument; in the event of a malfunction, discontinue use.

3. Application

3.1. Preparation

The medical devices must be tested prior to use with regard to the following aspects:

- Expiration Date.
- Undamaged packaging.
- Damage to the product (cracks on the tube, kinks, deformations).
- Check products for correct surfaces, correct assembly and functionality.
- Products that have not passed the functional test must not be used as their sterility and product safety are not guaranteed. Dispose of the products to be complained about accordingly or return them to the manufacturer.

3.1.1 Work Preparation & Testing

NOTE

Do not use products that have expired.

NOTE

If the packaging is damaged before use, replace the instrument with a new one.

NOTE

If you notice any irregularities, replace the instrument with a new one.

**WARNING**

Using a defective instrument can cause damage to the endoscope and result in injury to patients or users.

3.1.2 Compatibility and accessories

- Commercially available endoscopes of all types

3.2. Application

The brushes are carried out in uniformly axial dishes under permanent flushing with the cleaning agent permitted for the validated preparation of the endoscopes. The manufacturer's instructions must always be observed.

3.3. Preparation

**WARNING**

Reconditioning and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

4. Technical Data

The list of all variants as well as the materials and individual parts used can be found in the technical product data sheet and can be obtained on request from ENDO-FLEX GmbH the and/or on the website of the . ENDO-FLEX GmbH

5. Storage & Transportation

This product may only be transported and stored in the packaging provided for it. In addition, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store instruments near chemicals, disinfectants, or radioactive radiation.

SAFETY INSTRUCTIONS

Do not place objects on the storage packaging and the sterile barrier system!

6. Disposal

After use, this product may present a biological hazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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