Ballon-Dilatatoren

DE

EN

FR

IT

ES

NL

CZ

PT

SE

NO

DK

RU

CN

KR

TU

PL

SA

EL

HR

RO

DE	Ballon-Dilatator	Gebrauchsanweisung	3
EN	Dilation Balloon	Instruction Manual	11
FR	Ballon de dilatation	Manuel d'utilisation	18
IT	Dilatatore per palloncini	Manuale di istruzioni	25
ES	Dilatador de balón	Manual de instrucciones	32
NL	Ballon dilator	Handleiding	40
CZ	Balón dilatátor	Návod k použití	47
PT	Dilatador de balão	Manual de Instruções	54
SE	Ballong dilator	Instruktionsbok	61
NO	Ballong dilator	Instruksjonshåndbok	68
DK	Ballon dilator	Betjeningsvejledning	75
RU	диллятор баллона	Руководство по эксплуатации	82
CN	气球扩张器	说明书	90
KR	풍선 확장기	사용 설명서	97
TU	Balon dilatör	Kullanım talimatları	104
PL	Dilatatory balonem	Instrukcje użytkowania	111
SA	نولابلا عسوم	إدخستال تاميلعت	119
EL	Διαστολέας μπαλονιού	Οδηγίες χρήσης	126
HR	Balon dilator	Uputa za uporabu	134
RO	Dilatatoare cu balon	Instrucțiuni de utilizare	141

Symbole

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Symbols

	Medical Device		Working Length
	Catalog Number		Outside Diameter
	Batch		Minimum Working Channel
	Warning! Read the Instructions for use before use		Balloon Length / Balloon Diameter
	Sterilized with ethylene oxide		Maximum guide wire diameter
	Usable until		Only for the attending physician
	For single use only. Do not reuse		
	Do not re-sterilize!		
	Keep dry		
	Protect from sunlight		
	Manufacturer		
	Date Of Manufacture		
	Do not use if packaging is damaged		

Copyright

The copyright to these instructions for use remains with the company ENDO-FLEX GmbH. All rights, including those of translation, are reserved. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), in whole or in part, without the written permission of the company ENDO-FLEX GmbH or processed, reproduced, translated or distributed using electronic systems. Violations can result in criminal penalties.

Limitation Of Liability

The does ENDO-FLEX GmbH not assume any liability for personal injury, damage to property, damage to the product or consequential damage caused by non-compliance with these instructions for use, improper use of the product or any other action by unqualified and not ENDO-FLEX GmbH trained personnel. It is also forbidden to carry out modifications or technical modifications to the product by yourself.

Purpose of the Instructions for use

These Instructions for use are intended for use by the medical personnel working with the product. It describes the structure, function and operation of the product. Read these Instructions for use carefully and make sure that you understand them and that they are accessible to healthcare professionals

do it!

Follow the instructions to ensure safe working. Also, follow the instructions in other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The use of the product may only be carried out by qualified and, if necessary, on request ENDO-FLEX GmbH by trained medical personnel. The use of this product by healthcare professionals is always carried out taking into account the individual benefits of the patient at their risk.

Always have a copy of these Instructions for use at hand when using.

Classification

These Dilation Balloon are classified as Class **Is** products according to Medical Device Regulation (MDR) – (2017/745) via medical devices.

Reporting Obligation

The indicates ENDO-FLEX GmbH that the user is obliged to report to the competent authority of his country and at the same time all adverse events occurring during the use of this product and concerning this ENDO-FLEX GmbH product.

Language

These and all other instructions for use are available as a PDF download at www.endo-flex.de .

1. Product Description

The balloon dilators offered by ENDO-Flex GmbH generally have a distal transparent plastic balloon, a double-lumen plastic tube and two proximal luer-lock connectors, one of which serves as contrast medium and guide wire access and the second is provided with a 2-way tap. To fill the balloon with the filling medium and then block the filling volume of the balloon during treatment.


Due to the application, a distinction is made between biliary balloon dilators and achalasiathem balloon dilators.

Biliary Balloon Dilators:

These are placed via endoscope (TTS) by means of guide wire. The total length is 200 to 300cm with a tube diameter of CHF 7, a balloon diameter between 4 and 15mm and a balloon length between 20 and 30mm.

Achalasiathem Balloon Dilators:

These are placed without an endoscope (OTW) by means of a guide wire. The total length is 100cm for a tube diameter of CHF 14, a balloon diameter between 30 and 40mm and a balloon length of 80mm.

1	Balloon	2	Tube
3	Guide Wire Channel	4	Inflation channel with 2-way tap

This Instructions for use is valid for the products/product groups listed below:

Achalasia:

- 3413080
- 3413580
- 3414080

Biliary:

- 3410420
- 3410620
- 3410830
- 3411030
- 3411230
- 3410420-300
- 3410620-300
- 3410830-300

1.1. Contents of the packaging

- 1 Carton
- 1 Balloon Dilator SU (Single use) individually sterile packed
- 1 Instructions For Use

1.2. Intended Use

Balloon dilators are used for endoscopic dilatation of stenosis in the digestive tract. The application is done with or without an endoscope.

1.2.1 Duration Of Use

The application duration is specified for less than 60 minutes.

1.2.2 Lifetime

The life of the product is 3 years.

1.2.3 Indication

- Clinically relevant stenosis in esophagus, ductus choledochus
- Dysphagia of esophagus

1.2.4 Contraindications

- Non-sober patient
- Peritonitis, acute abdomen, e.g. Intestinal Perforation, Ileus
- Sepsis
- Comorbidity, such as severe cardiopulmonary disease and decompensation

- Uncontrollable hemorrhagic diathesis
- Pregnancy
- Recently created gastrointestinal anastomosis

1.2.5 Potential Complications/Risks

There is an increased risk of hypoxemia, hypercapnia, blood pressure drop, arrhythmias, and aspiration due to reduced protective reflexes when sedation is performed during endoscopic examination. Hypoxemia is also produced without sedation during endoscopic examinations due to the advance of the endoscope.

Possible injuries associated with endoscopic examinations may include: Perforations, bleeding, infections such as acute pancreatitis.

1.2.6 Target group of users

Trained specialists and specialists.

1.3. Patient Population

The patient or patient target group is derived from the indication of the responsible physician, who treats the patient in a diagnostic or therapeutic manner within the scope of an endoscopy (the leading procedure in itself), according to the intended use of the medical device. Limitations to patient population or target group are not known.

1.3.1 Use in minors

The user makes the use of the product in minors dependent on whether the physiological and anatomical conditions of the patient allow the use of the product.

1.3.2 Use in pregnant or breastfeeding women

The indication for the use of the product in pregnant or nursing women must be limited by the user on the basis of the individual physiological and anatomical conditions.

1.4. Sterility

These disposable medical devices are supplied sterile. Treatment and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

2. Safety Information

The safety and warning information used in these Instructions for use is based on ANSI Z535 standards.

DANGER	Warns of an imminent danger that leads to death or serious injury.
WARNING	Warns of a potentially dangerous situation that could lead to death or serious injury.
CAUTION	Warns of a potentially dangerous situation that could lead to moderate or minor personal injury.
NOTE	Warns of a potentially dangerous situation that could lead to property damage and environmental damage.
SAFETY INSTRUCTIONS	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

2.1. General

Follow these instructions, as well as instructions from the compatible components and hospital-internal regulations for infection prevention, safe use, cleaning and sterilization.

- Only sterile if the packaging is undamaged or unopened!
- For single use only! Do not reuse, reprocess, or sterilize multiple times. Re-use, reprocessing or repeated sterilization of the instrument may impair its structural integrity and may cause malfunctions, resulting in contamination, infections and serious injury.
- If the instrument accidentally becomes dirty before treatment, discard it immediately! Do not apply any cleaning agents!
- Do not use after expiry of the expiry date!
- All components should be carefully checked for compatibility and integrity prior to use. Do not use defective instruments! If any defects occur, discard the instrument and replace it with a new one.
- Never use the product outside the recommended technical specifications (intended use).
- Never manipulate the structural conditions of the instrument, avoid kinks and other damage, immediately discontinue use in the event of a malfunction!
- It is absolutely necessary to wear protective clothing (gloves, mouthguards, goggles, gown, etc.)!
- Compare the technical specifications of the product with those of the endoscope used. The working channel diameter must be at least 0.2 mm larger than the outer diameter of the instrument.
- Never force instruments into the working channel!

SAFETY INSTRUCTIONS

Take appropriate measures to protect the patient and yourself from contamination and injury. Follow national and local hygiene and safety guidelines. If in doubt, contact your local representative.

WARNING

Do not use defective instruments! If defects occur, discard the instrument and replace it with a new one.

NOTE

Never force instruments into the working channel.

NOTE

Never use the product outside the recommended technical specifications.

WARNING

Never manipulate the structural conditions of the instrument; in the event of a malfunction, discontinue use.

2.2. Specific

WARNING

Warning about general risks and complications with ERCP

- Post-ERCP bleeding
- Post-ERCP pancreatitis (PEP)
- Cholangitis and cholecystitis
- Perforation
- Spleen
- Long-term effects of sphincterotomy
- Air Embolism
- Abnormal response to sedative
- Lung and cardiovascular problems
- Tissue damage caused by X-ray radiation (especially in pediatric use)
- Infection due to inadequate treatment of endoscopes, instruments and accessories.

3. Application

3.1. Preparation

The medical devices must be tested prior to use with regard to the following aspects:

- Expiration Date.
- Undamaged packaging.
- Damage to the product (cracks on the tube, kinks, deformations).
- Check products for correct surfaces, correct assembly and functionality.
- Products that have not passed the functional test must not be used as their sterility and product safety are not guaranteed. Dispose of the products to be complained about accordingly or return them to the manufacturer.

3.1.1 Work preparation & Test

NOTE

If the packaging is damaged before use, replace the instrument with a new one.

NOTE

Do not use products that have expired.

NOTE

If you notice any irregularities, replace the instrument with a new one.

WARNING

Using a defective instrument can cause damage to the endoscope and result in injury to the patient or user.

1. Have a suitable inflation device with pressure gauge ready.

NOTE

Verify that the inflation device provides the necessary pressure and volume for inflation.

2. Provide adequate filling fluid (for example, a thin contrast agent, sterile, isotonic saline, or a mixture of both).
3. If you want to work with guide wire support, a suitable guide wire must be provided.
4. Open the package and carefully remove the catheter.
5. Check the catheter for damage.
6. Apply a slight vacuum to the Luer connector labeled "BALLOON" with a syringe and remove the protec-

tive cover.

- To improve the sliding properties in the endoscope, the balloon should be sprayed evenly with silicone spray.

NOTE

Only silicone should be used as a lubricant.

Test/Fill

To remove residual air from the balloon segment, the following procedure should be followed:

- Fill a 20ml syringe or other inflation device with approximately 5 ml of the selected fluid.
- After attaching the syringe or inflation device to the dilatation balloon filling lumen, orient the catheter vertically with the distal tip of the dilatation balloon pointing downwards.
- Apply vacuum and aspirate for at least 30 seconds. Slowly release the pressure and allow the contrast agent to fill the system.
- Switch to a fully filled syringe or a fully filled inflation device.
- Repeat aspiration without introducing air. Before inserting the dilation balloon into the endoscope, apply a vacuum in the inflator/balloon and maintain it when inserting it into the working channel.

NOTE

Significant amounts of residual air can lead to uneven expansion of the balloon. If this is noticed, repeat the filling steps.

3.1.3 Compatibility and accessories

- Endoscopes with a working channel from 2.8 mm Ø
- Inflation device with Luer-Lock connector and pressure display
- Guide wires up to .035"

3.2. Application

Biliary Dilation:

- Push the balloon dilator into the endoscope with short (2 to 3 cm) movements

NOTE

Due to the different endoscope versions, a slight resistance can be felt immediately upon entry into the endoscope and about 2 to 3 cm before exit at the distal end of the working channel.

- Push the balloon under endoscope control to the desired position.
- Filling the balloon with a suitable inflation device under constant pressure control. Sterile, isotonic saline, thin-fluid X-ray contrast agents or mixtures

of both are suitable as medium.

WARNING

Never fill balloon with air or other gaseous media!

NOTE

Do not exceed the specified maximum inflation pressure. If the balloon bursts or a loss of pressure occurs, the balloon must be completely emptied and carefully removed together with the endoscope. Do not attempt to pull a burst and/or broken balloon through the endoscope. Continue the procedure with a new dilation balloon.

- After reaching the desired diameter, adjust the balloon pressure (do not exceed maximum inflation pressure) as required.

NOTE

After a working pressure has been set at the inflator unit, the subsequent pressure distribution in the catheter system results in a small pressure drop at the pressure display, which should be stable after a few seconds.

- The pressure is maintained until the desired dilation is achieved.

NOTE

During dilation, the pressure gauge may fluctuate.

- After application, empty the balloon. This takes about 20 to 30 seconds depending on the balloon size and filling medium.
- For complete emptying, the proximal end of the balloon must remain endoscopic while the filling device is used to apply vacuum.

NOTE

Do not retract the catheter until the balloon is completely empty

- By straightening the distal endoscope end, the removal of the catheter is simplified. Any superfluous curvature in the working channel makes it difficult to pull out the catheter through the endoscope.
- Slowly pull the dilation balloon out through the endoscope.

NOTE

If excessive resistance is felt, remove the endoscope and balloon assembly to avoid tissue injury or endoscope damage.

- After use, the instrument must be disposed of in accordance with the relevant legal and internal clinical

regulations.

Agalasic dilation:

1. Push the balloon to the desired position.
2. Filling the balloon with a suitable inflation device under constant pressure control. Sterile, isotonic saline, thin-fluid X-ray contrast agents or mixtures of both are suitable as medium.

⚠ WARNING

Never fill balloon with air or other gaseous media!

NOTE

Do not exceed the specified maximum inflation pressure. If the balloon bursts or a loss of pressure occurs, the balloon must be completely emptied and carefully removed. Continue the procedure with a new dilation balloon.

3. After reaching the desired diameter, adjust the balloon pressure (do not exceed maximum inflation pressure) as required.

NOTE

After a working pressure has been set at the inflator unit, the subsequent pressure distribution in the catheter system results in a small pressure drop at the pressure display, which should be stable after a few seconds.

4. The pressure is maintained until the desired dilation is achieved.

NOTE

During dilation, the pressure gauge may fluctuate.

5. After application, empty the balloon. This takes about 20 to 30 seconds depending on the balloon size and filling medium.

NOTE

Do not retract the catheter until the balloon is completely empty.

After use, the instrument must be disposed of in accordance with the relevant legal and internal clinical regulations.

3.3. Preparation

⚠ WARNING

Reconditioning and sterilization is no longer possible after use and is prohibited! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

4. Technical Data

The list of all variants as well as the materials and individual parts used can be found in the technical product data sheet and can be obtained on request from ENDO-FLEX GmbH and/or on the website of the ENDO-FLEX GmbH

5. Storage & Transportation

This product may only be transported and stored in the packaging provided for it. In addition, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store instruments near chemicals, disinfectants, or radioactive radiation.

SAFETY INSTRUCTIONS

Do not place objects on the storage packaging and the sterile barrier system!

6. Disposal

After use, this product may present a biological hazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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