# Saug-Polypenfalle SU

DE

EN

FR

IT

ES

NL

CZ

PT

SE

NO

DK

RU

CN

KR

TU

PL

SA

EL

HR

RO

<b>DE</b>	Saug-Polypenfalle SU (Einmalgebrauch)	Gebrauchsanweisung	<b>3</b>
<b>EN</b>	Suction polyp trap SU	Instruction Manual	<b>8</b>
<b>FR</b>	Piège à polype par aspiration SU	Manuel d'utilisation	<b>13</b>
<b>IT</b>	Trappola per polipi ad aspirazione SU	Manuale di istruzioni	<b>18</b>
<b>ES</b>	Trampa de pólipos por succión SU	Manual de instrucciones	<b>23</b>
<b>NL</b>	Zuigpoliepenvanger SU	Handleiding	<b>28</b>
<b>CZ</b>	Odsávací polygová past SU	Návod k použití	<b>33</b>
<b>PT</b>	Suction polyp trap SU	Manual de Instruções	<b>38</b>
<b>SE</b>	Sugpolypfälla SU	Instruktionsbok	<b>43</b>
<b>NO</b>	Sugekopp i polypp SU	Instruksjonshåndbok	<b>48</b>
<b>DK</b>	Sugepolypfælde SU	Betjeningsvejledning	<b>53</b>
<b>RU</b>	Всасывающая ловушка для полипов SU	Руководство по эксплуатации	<b>58</b>
<b>CN</b>	抽吸式息肉捕集器 SU	说明书	<b>63</b>
<b>KR</b>	흡입 폴립 트랩 SU	사용 설명서	<b>68</b>
<b>TU</b>	Emme polip tutucu su	Kullanım talimatları	<b>73</b>
<b>PL</b>	Ssący separator polipów SU	Instrukcje użytkowania	<b>78</b>
<b>SA</b>	طفششلا دصم SU	مدخست الال تامي لعت	<b>83</b>
<b>EL</b>	Παγίδα πολυπόδων αναρρόφησης SU	Οδηγίες χρήσης	<b>88</b>
<b>HR</b>	Usisni polip zamka SU	Uputa za uporabu	<b>93</b>
<b>RO</b>	Capcană de aspirație pentru polipi SU	Instrucțiuni de utilizare	<b>98</b>

## Symbole

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Catalog number

Charge

Attention! Read the instructions for use before use

Usable until

For single use only. Do not reuse

Do not sterilize again!

Keep dry

Protect from sunlight

Manufacturer

Date of manufacture

Do not use if packaging is damaged

Only for the attending physician

Non-sterile

## Copyright

The copyright to these operating instructions remains with the company ENDO-FLEX GmbH. All rights reserved, including those of translation. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), either in whole or in part, without the written permission of the company ENDO-FLEX GmbH reproduced or processed, duplicated, translated or distributed using electronic systems. Violations may result in criminal penalties.

## Limitation of liability

The ENDO-FLEX GmbH assumes no liability for personal injury, property damage, damage to the product or consequential damage resulting from failure to observe these operating instructions, from improper use of the product or from the actions of unqualified personnel or personnel not trained by the ENDO-FLEX GmbH trained personnel. It is also prohibited to make unauthorized modifications or technical changes to the product.

## Purpose of the instructions for use

These instructions for use are intended for medical personnel working with the product. It describes the structure, function and operation of the product. Read these instructions for use carefully and make sure you understand them and make them available to the healthcare professional

make!

Follow the instructions to ensure safe working. Also observe the information in the instructions of other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The product may only be used by qualified medical personnel, who may be instructed by the ENDO-FLEX GmbH instructed medical personnel. The use of this product by healthcare professionals always involves weighing the individual patient's benefit against the patient's risk.

Always keep a copy of these instructions for use handy during use.

## Classification

The Suction polyp trap SU is classified as class I -Product according to Medical Device Regulation (MDR) - (2017/745) on medical devices.

## Obligation to report

The ENDO-FLEX GmbH points out that the user is obliged to report all adverse incidents occurring during the use of this product and concerning it to the competent authority of his country and at the same time to ENDO-FLEX GmbH report.


## Language

These and all other instructions for use are available at [www.endo-flex.de](http://www.endo-flex.de) as a PDF download.

## 1. Product description

The Suction polyp trap SU consists of four separate chambers for the collection of polyps and four open areas for direct aspiration.

The tightly closing lid is equipped with two nozzles, on each of which a suction hose can be attached.

- |   |  |
|---|--|
| 1 | Suction hose (left: to suction system & right: to endoscope) |
| 2 | Collecting and transport container with 4 chambers           |

These instructions for use are valid for the products listed below:

- **Suction polyp trap SU:** 900334H

### 1.1. Packing contents

The Suction polyp trap SU are supplied in packaging units of 6 pieces:

- 1 outer carton
- 6 Suction polyp trap SU individually unsterile packed
- 1 Instruction manual

### 1.2. Intended use

The Suction polyp trap SU serves as a collection and transport container for polyps removed during endoscopy for histological examination.

#### 1.2.1 Application duration

The application time is specified for less than 60 minutes.

#### 1.2.2 Lifetime

The service life of the product is 5 Years.

#### 1.2.3 Indication

- Endoscopic removal of ablated polyps

### 1.2.4 Contraindications

#### Contraindication

- There are no known contraindications for this product.

### 1.2.5 Potential complications/risks

- Under sedation during endoscopic examination, there is an increased risk of hypoxemia, hypercapnia, hypotension, arrhythmias, and aspiration due to decreased protective reflexes. Hypoxemia occurs even without sedation during endoscopic examinations due to the advancement of the endoscope.
- Possible injuries associated with endoscopic examinations may include: Perforations, bleeding, infections such as acute pancreatitis.

### 1.2.6 User target group

Trained specialists and professional staff.

### 1.3. Patient population

The patient or patient target group is derived from the indication of the responsible physician, who treats the patient diagnostically or therapeutically within the scope of an endoscopy (the leading intervention per se), in accordance with the intended use of the medical device. Limitations to patient populations or patient target groups are not known.

### 1.4. Sterility

These single-use medical devices are supplied

non-sterile. Reprocessing and sterilization is no longer possible and prohibited after use! The product may only be used once and must be disposed of after use, in accordance with national and internal guidelines.

## 2. Safety instructions

The safety and warning information used in these operating instructions is based on ANSI Z535 standards.

<b>DANGER</b>	Warns of an imminent danger that will result in death or serious bodily injury.
<b>WARNING</b>	Warns of a potentially dangerous situation that could result in death or serious bodily injury.
<b>CAUTION</b>	Warns of a potentially dangerous situation that could result in moderate or minor bodily injury.
<b>NOTE</b>	Warns of a potentially dangerous situation that could lead to property and environmental damage.
<b>SAFETY INSTRUCTION</b>	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

### 2.1. General

#### SAFETY INSTRUCTIONS

Protect the patient and yourself from contamination and injury by taking appropriate measures. Follow national and internal hygiene and safety guidelines. If in doubt, contact your local responsible person.

#### WARNING

Do not use defective instruments! If defects occur, discard the instrument and replace it with a new one.

#### NOTE

Never use the product outside the recommended technical specifications.

#### WARNING

Never tamper with the structural features of the instrument; abort the operation in the event of a malfunction.

## 2.2. Specific

### SAFETY INSTRUCTIONS

Do not attempt to remove the pre-assembled hose from the trap. This could lead to contact with the contents of the trap during removal.

## 3. Application

### 3.1. Work Preparation & Test

#### WARNING

Using a defective instrument can cause damage to the endoscope and injury to the patient or user.

#### NOTE

If the packaging is already damaged before use, replace the instrument with a new one.

#### NOTE

Do not use damaged packaging or products with expired expiration dates.

#### NOTE

If you notice any irregularities, replace the instrument with a new one.

### 3.1.1 Compatibility and accessories

Endoscopes / suction systems

## 3.2. Application

1. Check the instrument for safe operation and irregularities before the procedure.
2. If you notice any irregularities, replace the instrument with a new one.

### 3.2.1 Components

The suction polyp trap consists of:

- Vessel with four filter traps (labeled 1, 2, 3, and 4) for picking up polyps and four open areas (labeled 0) for direct aspiration.
- Cover with one hose each (including connection to the suction system) and one hose for connection to the endoscope.

### 3.2.2 Application

1. Connect the hose from the middle connection of the trap to the suction system (Fig. 1).
2. Connect the outer connector to the endoscope (Fig. 1).


Abb. 1

- For normal suction, turn the vessel until the arrow on the lid points to an open position ("0") (Fig. 2)


Abb. 2

- To pick up the first polyp, rotate the vessel until the arrow points to the number "1" (Fig. 3).


Abb. 3

- After picking up the polyp, repeat step 3 as needed for normal aspiration.
- Repeat step 4 for positions 2, 3 and 4 for each additional polyp (only for the same patient).
- Remove the hoses from the endoscope and suction system.
- Fill the suction polyp trap through the center port or after removing the lid with a fixative. After filling with fixative, close the polyp trap again by pressing the lid firmly.
- Connect both hose ends firmly with the hose connector. The suction polyp trap is thus sealed (Fig. 4).


Abb. 4

### 3.3. Preparation

**WARNING**

Reprocessing and sterilization is no longer possible and prohibited after use! The product may only be used once and must be disposed of after use in accordance with national and local guidelines

### 4. Technical data

The list of all variants as well as the materials and individual parts used can be found in the Technical Product Data Sheet and is available on request from the ENDO-FLEX GmbH and/or on the website of ENDO-FLEX GmbH available.

### 5. Storage & Transport

This product may only be transported and stored in the packaging provided. Beyond that, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store the instruments near chemicals, disinfectants or radioactive radiation.

**SAFETY INSTRUCTIONS**

Do not place any objects on the storage packaging and the sterile barrier system!

### 6. Disposal

After use, this product may present a biohazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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