

EasyCollect – Polyp & Foreign Body Retrieval Devices



Oval Thermo-Formed Plastic Bag

Latex-free

Rotatable

360°

Bag with Perforated Reservoir

For drainage of liquids

EasyCollect – Polyp & Foreign Body Retrieval Devices – Single Use



EasyCollect – Polyp and Foreign Body Retrieval Device

Oval thermo-formed plastic bag

Catalogue No.	Bag Width (mm)	Bag Depth (mm)	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE1518-B15	15	12	1.8	160	2.0
NE1522-G25	25	17	2.6	230	2.8
NE1522-G35	35	20	2.6	230	2.8

LEF-GFB-SU-00-EN / LL-0361 Vers. 1.0 (07/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de