

Biliary Dilation Catheters

Single Use



Radiopaque Teflon Tube with Distal Metal Ring

Precise position control under X-ray vision

Lateral Luer-Lock Connector for Injection of Contrast Medium

Regular 30 mm long Tapered Tip

Controlled and safe dilatation of the stenosis

Biliary Dilation Catheters – Single Use



Biliary Dilation Catheter

Working Length 200 cm; radiopaque teflon tube with distal metal ring; compatible with guide wire

Catalogue No.	Diameter Tubus (French)	Diameter Tip (French)	Compatible with Guide Wire (inch)	Min. Accessory Channel
E20000053	5	3	.025	2.0
E20000064	6	4	.035	2.2
E20000075	7	5	.035	2.8
E20000085	8	5	.035	3.2
E20000095	9	5	.035	3.2
E20000105	10	5	.035	3.7
E20000107	10	7	.035	3.7
E20000127	12	7	.035	4.2

LEF-GDC-SU-00-EN / LL-0393 Vers. 1.0 (11/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de