

# Lithotripsy Baskets and Spirals

Single Use



LEF-GML-SU-00-EN / LL-0385 Vers. 4.0 (04/2021)

## Combination Concept

Components for planned and emergency lithotripsies

## Simple and Safe Handling

Convenient handling for the assembly and detachability of components

## Pre-assembled Set Configurations

For specific treatment requirements



### ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany  
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11  
Email: [sales@endo-flex.de](mailto:sales@endo-flex.de) • [www.endo-flex.de](http://www.endo-flex.de)

# Lithotripsy Baskets – Single Use


# Lithotripsy Baskets – Single Use



## Lithotripsy Basket – 4 Braided Wire

Suitable for lithotripsy handle type I & III; tube length: 195 cm; total system length: 400 cm

Catalogue No.	Description	Sheath Outer Diameter (mm)	Suitable for Lithotripsy Spiral Ø (mm)	Basket Width in Opened Condition (mm)	Basket Length in Stretched Condition (mm)
E122264N	4-Wire Basket Type	2.6	2.6	~20	40
E122265N	4-Wire Basket Type	2.6	2.6	~25	50
E122266N	4-Wire Basket Type	2.6	2.6	~30	60
E122267N	4-Wire Basket Type	2.6	2.6	~35	70
E122336N	4-Wire Basket Type	2.9	3.3	~30	60



## Lithotripsy Basket – 6 Braided Wire


Suitable for lithotripsy handle type I & III; tube length: 195 cm; total system length: 400 cm

Catalogue No.	Description	Sheath Outer Diameter (mm)	Suitable for Lithotripsy Spiral Ø (mm)	Basket Width in Opened Condition (mm)	Basket Length in Stretched Condition (mm)
E124263N	6-Wire Dormia Type	2.6	2.6	~12	30
E124264N	6-Wire Dormia Type	2.6	2.6	~18	40
E124265N	6-Wire Dormia Type	2.6	2.6	~25	50
E124266N	6-Wire Dormia Type	2.6	2.6	~30	60
E124267N	6-Wire Dormia Type	2.6	2.6	~35	70
E124335N	6 Wire Dormia Type	2.9	3.3	~25	50
E124337N	6-Wire Dormia Type	2.9	3.3	~35	70
E128336N	6-Wire Basket Type	2.9	3.3	~30	60



# Lithotripsy Spirals – Single Use

# Lithotripsy Handles – Reusable



## Lithotripsy Spiral

Made of stainless steel

Catalogue No.	Suitable for Lithotripsy Handle	Outer Diameter (mm)	Total Length (cm)
E120261	Type I	2.6	195
E120602	Type III	2.6	195



## Emergency Lithotripsy Spiral

Made of stainless steel

Catalogue No.	Suitable for Lithotripsy Handle	Outer Diameter (mm)	Total Length (cm)
E120604	Type III	4.0	100




## Lithotripsy Set Type III

Consists of storage case, handle, spirals and baskets

Catalogue No.	Suitable for Lithotripsy Handle	Outer Diameter (mm)	Component Catalogue No.
E120610	Type III	1 x Lithotripsy Handle Type III	120600
		1 x Lithotripsy Storage Case	120601
		3 x Lithotripsy Spirals	E120602
		1 x Emergency Lithotripsy Spiral	E120604
		1 x Lithotripsy Basket	E122265N
		1 x Lithotripsy Basket	E122266N
		1 x Lithotripsy Basket	E122267N

## Fitting Lithotripsy Handles – Mechanical Lithotripsy

Catalogue No.	Description
120000N	<p><b>Lithotripsy Handle Type I</b></p> <p>Reusable, with coaxial fixation, ergonomic hand wheel suitable for planned Lithotripsy only</p> 
120600	<p><b>Lithotripsy Handle Type III</b></p> <p>Reusable, user-friendly twist &amp; pull technology, 1-push safety system suitable for planned and emergency Lithotripsy</p> 