

# Retrieval Net

Standard & Rotatable Retrieval Net



**Used for endoscopic removal of stones  
and foreign bodies in the GI tract.**

# Retrieval Net

## Standard Retrieval Net

Catalog No.	Sales Unit	Outer Ø	Working Length	Requested min. Channel	Net Ø
EFB-11-25-240	Box of 10	2.6mm	240cm	2.8mm	25mm
EFB-11-32-240	Box of 10	2.6mm	240cm	2.8mm	32mm
EFB-11-25-180	Box of 10	2.6mm	180cm	2.8mm	25mm
EFB-11-32-180	Box of 10	2.6mm	180cm	2.8mm	32mm

## Rotatable Retrieval Net

Catalog No.	Sales Unit	Outer Ø	Working Length	Requested min. Channel	Net Ø
EFB-12-25-240	Box of 10	2.6mm	240cm	2.8mm	25mm
EFB-12-32-240	Box of 10	2.6mm	240cm	2.8mm	32mm
EFB-12-25-180	Box of 10	2.6mm	180cm	2.8mm	25mm
EFB-12-32-180	Box of 10	2.6mm	180cm	2.8mm	32mm



- Function of rotation is designed to accommodate physician technique preference.
- Nylon material engineered for better supporting strength.
- Net designed to obtain a sizeable amount of tissue.

### ENDO-FLEX GmbH

Alte-Hünxer-Str. 115  
46562 Voerde, Germany  
Telefon: +49 281 94400-0  
Fax: +49 281 94400-11  
Email: [sales@endo-flex.de](mailto:sales@endo-flex.de)  
[www.endo-flex.de](http://www.endo-flex.de)