

# Bite Blocks

Single Use – Reusable

## Bite Blocks

- Ergonomic Adaption to Patient's Facial Outline  
*High Patient Comfort*
- Complete Latex-Free   
*Eliminates Risk of Latex Allergy Reactions*
- Stackable Dispenser  
*Hygienic and User-Friendly Storage*



# Bite Blocks

Single Use – Reusable

## Single Use

### Bite Block Types



GAC-01-01-001



GAC-01-03-001



GAC-01-03-002

### Bite Block – Complete Latex-Free, Individually Packed in a Dispenser Box

Catalog No.	Description	Opening Size	Unit
GAC-01-01-001	With Tongue Pressure and Rubber Headstrap	W 19 mm x H 18 mm	Box of 50
GAC-01-03-001	With Rubber Head Strap	W 26 mm x H 18 mm	Box of 50
GAC-01-03-002	With Textile Head Strap	W 23 mm x H 20 mm	Box of 100



### Bite Block Type

### Bite Block – Complete Latex-Free, Individually Packed in a Dispenser Box

Catalog No.	Description	Opening Size	Unit
GAC-01-01-002	Without Head Strap	W 26 mm x H 18 mm	Box of 100



## Reusable

### Bite Block Type

### Bite Block – Complete Latex-Free, Individually Packed

Catalog No.	Description	Opening Size	Unit
GAC-01-02-000	Without Head Strap	W 26 mm x H 18 mm	Box of 5