


# Cytology Brushes

Single Use

## Cytology Brushes

- Brush with Metal Tip  
*Efficient tube closure to secure smear yield*
- Set Versions including Brush and Scissor  
*Safe and user-friendly specimen handling*
- Sophexa Easy Brush  
*Guide wire compatible cytology brush for smear acquisition from the bile duct*



# Cytology Brushes

Single Use

## Bronchoscopy

### Cytology Brush – with Conical Metal Tip and Ergonomic Handle – sterile packed

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-120	Cytology Brush	120 cm	1.8 mm	2.0 mm
GCB-01-24-120	Cytology Brush	120 cm	2.45 mm	2.8 mm
GCB-02-18-120	Set incl. Cytology Brush & Scissor	120 cm	1.8 mm	2.0 mm

## Gastroscopy

### Cytology Brush – with Conical Metal Tip and Ergonomic Handle – sterile packed

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-180	Cytology Brush	180 cm	1.8 mm	2.0 mm
GCB-01-24-180	Cytology Brush	180 cm	2.45 mm	2.8 mm
GCB-02-18-180	Set incl. Cytology Brush & Scissor	180 cm	1.8 mm	2.0 mm

## Colonoscopy

### Cytology Brush – with Conical Metal Tip and Ergonomic Handle – sterile packed

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-01-18-230	Cytology Brush	230 cm	1.8 mm	2.0 mm
GCB-01-24-230	Cytology Brush	230 cm	2.45 mm	2.8 mm

## Sophexa Easy Brush

Guide Wire Compatible Cytology Brush for Smear Acquisition mainly from the Bile Duct.

### Cytology Brush – with Ellipsoid Tip and Ergonomic Handle – sterile packed

Catalog No.	Description	Working Length	Sheath Outer Diameter	Minimum Accessory Channel Diameter
GCB-03-24-230	Cytology Brush	230 cm	2.45 mm	2.8 mm