

Nasal Biliary Drainages


Single Use

Nasal Biliary Drainages

- Various Types and Sizes Available
Fit for all anatomic cases
- Well Balanced Wall Diameter
Maximum inner lumen with highest form stability & kink resistance
- Tapered & Round Ends
Atraumatic implantation




Nasal Biliary Drainages



Nasal Biliary Drainages Sets

Components: Drainage Catheter with Position Marker, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Connecting Tube
		French Size	Length	Length	Length
GND-02-05-290S	Standard	5.0	250 cm	50 cm	135 cm
GND-02-06-290S	Standard	6.0	250 cm	50 cm	135 cm
GND-02-07-290S	Standard	7.0	250 cm	50 cm	135 cm
GND-02-10-290S	Standard	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter with Position Marker, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Connecting Tube
		French Size	Length	Length	Length
GND-03-05-290S	Standard with Pigtail	5.0	250 cm	50 cm	135 cm
GND-03-06-290S	Standard with Pigtail	6.0	250 cm	50 cm	135 cm
GND-03-07-290S	Standard with Pigtail	7.0	250 cm	50 cm	135 cm
GND-03-10-290S	Standard with Pigtail	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Connecting Tube
		French Size	Length	Length	Length
GND-01-05-290S	Pigtail	5.0	250 cm	50 cm	135 cm
GND-01-06-290S	Pigtail	6.0	250 cm	50 cm	135 cm
GND-01-07-290S	Pigtail	7.0	250 cm	50 cm	135 cm
GND-01-10-290S	Pigtail	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Connecting Tube
		French Size	Length	Length	Length
GND-04-05-290S	O – Type	5.0	250 cm	50 cm	135 cm
GND-04-06-290S	O – Type	6.0	250 cm	50 cm	135 cm
GND-04-07-290S	O – Type	7.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Connecting Tube
		French Size	Length	Length	Length
GND-05-05-290S	J – Type	5.0	250 cm	50 cm	135 cm
GND-05-06-290S	J – Type	6.0	250 cm	50 cm	135 cm
GND-05-07-290S	J – Type	7.0	250 cm	50 cm	135 cm


All Drainage Catheters are also available individually (Catalog set no. without “S”).

Nasal Biliary Drainages


Single Use

Individual Components

Nasal Tube

Catalog No.	French Size	Length	Product
GND-41-05-050	5.0	50 cm	
GND-41-06-050	6.0	50 cm	
GND-41-07-050	7.0	50 cm	
GND-41-10-050	10.0	50 cm	

Multi-Medfix Injection Connector

Catalog No.	French Size	Product
GND-43-05-000	5.0	
GND-43-06-000	6.0	
GND-43-07-000	7.0	
GND-43-10-050	10.0	

Connecting Tube

Catalog No.	Length	Product
GND-42-03-135	135 cm	