


# Stone Extraction Baskets

Single Use

## Stone Extraction Baskets

- 45° Injection Port & Integrated Sealing  
*User-friendly and easy injection of contrast medium*
- Detachable Instrument Holder & Outer Sheath  
*Applicable to emergency lithotripsy spiral and StoneBuster handle in emergency cases*
- Ergonomic Instrument Holder  
*Simple and safe basket manipulation*



# Stone Extraction Baskets

Single Use

## Features & Benefits

- ↗ 45° Injection Port & Integrated Sealing  
*User-friendly and easy injection of contrast medium*
- ↗ Detachable Instrument Holder & Outer Sheath  
*Applicable to emergency lithotripsy spiral and StoneBuster handle in emergency cases*
- ↗ Ergonomic Instrument Holder  
*Simple and safe basket manipulation*

### Stone Extraction Basket Set – Set Consists of Basket, Outer Sheath (Working Length 180 cm), Sealing and Instrument Holder

Catalog No.	Description	Basket Diameter	Basket Length	Outer Diameter
GEB-11-19-420	4-Wire Baskets	20 mm	42 mm	1.9 mm
GEB-11-19-425	4-Wire Baskets	25 mm	52 mm	1.9 mm
GEB-11-19-430	4-Wire Baskets	30 mm	61 mm	1.9 mm
GEB-11-19-620	6-Wire Baskets	20 mm	42 mm	1.9 mm
GEB-11-19-625	6-Wire Baskets	25 mm	52 mm	1.9 mm
GEB-11-19-630	6-Wire Baskets	30 mm	61 mm	1.9 mm
GEB-11-25-420	4-Wire Baskets	20 mm	42 mm	2.5 mm
GEB-11-25-425	4-Wire Baskets	25 mm	52 mm	2.5 mm
GEB-11-25-430	4-Wire Baskets	30 mm	61 mm	2.5 mm
GEB-11-25-620	6-Wire Baskets	20 mm	42 mm	2.5 mm
GEB-11-25-625	6-Wire Baskets	25 mm	52 mm	2.5 mm
GEB-11-25-630	6-Wire Baskets	30 mm	61 mm	2.5 mm

Distributed by:

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany  
☎ +49 80 32 / 973- 379 | 📠 +49 80 32 / 973- 339 | ✉ sales@medi-globe.de

[www.medi-globe.de](http://www.medi-globe.de)

L-GEB-SU-00-ENG Ver. 6.0 (2015-11)