


ERCP Catheters – Single Use



ERCP Catheter – with Standard Stylet

With colored markers; compatible with guide wire .035"

Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E4318-GW	Standard	1.8	200	2.0
E4418-GW	Tapered	1.8	200	2.0
E4518-GW	Conical Metal	1.8	200	2.0
E4618-GW	Tapered Metal	1.8	200	2.0
E4818-GW	Hour Glass Metal	1.8	200	2.0



Standard Tip

ERCP Catheter – with Standard Stylet

Compatible with guide wire .035"

Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E43162FD35	Standard	1.65	200	2.0



LEF-GER-SU-00-EN / LL-0371 Vers. 1.0 (10/2020)

ERCP Catheters

Single Use



Guide Wire Compatibility

Color-Marked

Versions with Metal Tip




ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de

ERCP Catheters – Single Use

ERCP Catheters – Single Use



ERCP Catheter – with Luer-Lock Stylet

With colored markers, compatible with guide wire .035"


Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E4318	Standard	1.8	200	2.0
E4418	Tapered	1.8	200	2.0



ERCP Catheter – with Standard Stylet

With colored markers, compatible with guide wire .035"

Catalogue No.	Tip Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E4318-GW-250	Standard	1.8	250	2.0
E4318-GW-280	Standard	1.8	280	2.0



ERCP Catheter – with Luer-Lock-Stylet

Contrast medium injection while guide wire is inserted; with colored markers, three access ports, compatible with guide wire .035"

Catalogue No.	Tip Shape	Catheter Outer Diameter Prox/ Dist (mm)	Working Length (cm)	Min. Accessory Channel (mm)
E4422-18FD	Tapered	2.45/1.8	200	2.8
E4422-18FD-280*	Tapered	2.45/1.8	280	2.8
E4522-18FD	Conical Metal	2.45/1.8	200	2.8
E4622-18FD-230	Tapered Metal	2.45/1.8	230	2.8
E4722-18FD	Round Metal	2.45/1.8	200	2.8
E4822-18FD	Hour Glass Metal	2.45/1.8	200	2.8



* Packing unit 1 pc.