

Nasal Biliary Drainages

Single Use



Various Types and Sizes

Taylored fit for all anatomic cases


Balanced Wall Diameter

Maximum inner lumen with form stability and kink resistance

Tapered and Round Endings

Atraumatic placement

Nasal Biliary Drainages – Single Use



Nasal Biliary Drainage

Radiopaque probe with position marking (probe length 250 cm); nasal catheter (length 50 cm); drainage connection tube with stop-cock and Luer-Lock connector

Catalogue No.	Catheter Shape	Probe Size (French)	Drainage Connection Tube Length (cm)	Compatible Guide Wire (inch)	Min. Accessory Channel (mm)
206185	Classic	5	135	.035	2.0
206225	Classic	7	135	.035	2.8
206305	Classic	10	135	.047	3.7
207185	Classic + Pigtail	5	135	.035	2.0
207225	Classic + Pigtail	7	135	.035	2.8
207305	Classic + Pigtail	10	135	.047	3.7
208185	Pigtail	5	135	.035	2.0
208225	Pigtail	7	135	.035	2.8
208305	Pigtail	10	135	.047	3.7

LEF-GND-SU-00-EN / LL-0337 Vers. 3.0 (10/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de